

Product Datasheet

Anti-Calreticulin Antibody (orb1148634)

Description

Anti-Calreticulin Antibody Rabbit Polyclonal

Species/Host	Rabbit
Reactivity	Bovine, Canine, Gallus, Guinea pig, Hamster, Human, Monkey, Mouse, Porcine, Rabbit, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human calreticulin synthetic peptide conjugated to KLH. Host Species: Rabbit Antibody Class: Polyclonal Preservatives: None.
Target	Calreticulin
Preservatives	Whole Antiserum
Form/Appearance	Liquid
Storage	This antibody is stable for at least one (1) year at -20C. Avoid repeated freezing and thawing.
Note	For research use only
Clonality	Polyclonal
Purity	Serum
Uniprot ID	P27797
Dilution Range	Immunoblotting: use at 1:5,000-1:10,000 dilution. A band of ~63 kDa is detected. These are recommended concentrations. User should determine optimal concentrations for their application. Positive control: HeLa cell lysate.
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558