

Product Datasheet

Anti-RPA32 Antibody (orb1148667)

Description	Anti-RPA32 Mouse Monoclonal Antibody
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	N-terminally 6x-His tagged fusion protein corresponding to full-length human RPA32 expressed in E. coli. Host Species: Mouse Antibody Class: IgG2a/ α Stabilizers: None Preservatives: None. Available on request.
Target	RPA32
Preservatives	PBS, pH 7.4
Form/Appearance	Liquid
Storage	This antibody is stable for at least one (1) year at -70C. Avoid multiple freeze- thaw cycles.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P15927
Dilution Range	Immunoblotting, immunoprecipitation: use at 0.2-1 ug/ml. Positive control: HeLa, Raji, and 3T3 cells and recombinant protein.
Expiration Date	12 months from date of receipt.