

Product Datasheet

Anti-Rad50 Antibody (orb1148670)

Description	Anti-Rad50 Mouse Monoclonal Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	GST fusion protein corresponding to the complete coding region (aa 1-425) of Rad50 expressed in E. coli. Host Species: Mouse Antibody Class: IgG1/ μ Stabilizers: None Preservatives: None. Available on request.
Target	Rad50
Preservatives	PBS, pH 7.4
Form/Appearance	Liquid
Storage	This antibody is stable for 1 month at 4°C or one year at -20°C. Avoid multiple freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Uniprot ID	Q92878
Dilution Range	Immunoblotting, immunoprecipitation: use at 0.1-2 ug/ml. Positive control: Raji cells.
Expiration Date	12 months from date of receipt.