

Product Datasheet

Anti-Mre11 Antibody (orb1148671)

Description	Anti-Mre11 Mouse Monoclonal Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	GST fusion protein corresponding to aa 182-582 of Mre11.
Target	Mre11
Preservatives	PBS, pH 7.4.
Form/Appearance	Liquid
Storage	This antibody is stable for at least one (1) year at -20C. Avoid multiple freeze- thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P49959
Dilution Range	Immunoblotting: Use at 0.5-2ug/ml. A band of ~80kDa is detected. Detection of Mre11 in 30ug of whole cell extracts. Immunofluorescence: use at 1-5ug/ml. Detection of Mre11 in HeLa cells fixed in 4% paraformaldehyde for 15 minutes. These are recommended concentrations. Endusers should determine optimal concentrations for their applications.
Expiration Date	12 months from date of receipt.