

Product Datasheet

Anti-Id3 Antibody (orb1148673)

Description	Anti-Id3 Mouse Monoclonal Antibody
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Mouse Id3 gst-fusion protein (aa 1-119) expressed in E. coli. Host Species: Mouse Stabilizers: None Preservatives: None. Available on request.
Target	Id3
Preservatives	10 mM PBS, pH 7.4.
Form/Appearance	Liquid
Storage	This antibody is stable for at least one (1) year at -20C. Avoid multiple freeze- thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Uniprot ID	Q02535
Dilution Range	Immunoblotting and Immuno- precipitation: use at 1-5 ug/ml. Positive controls: Id3 transfected Sol8 muscles cells and recombinant fusion protein.
Expiration Date	12 months from date of receipt.