

Product Datasheet

Anti-Fen-1 Antibody (orb1148674)

Description	Anti-Fen-1 Mouse Monoclonal Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant human FEN-1 protein (aa 1-380) purified from baculovirus-infected Sf-9 cells. Host Species: Mouse Stabilizers: None Preservatives: None. Available on request.
Target	Fen-1
Preservatives	10 mM PBS, pH 7.4.
Form/Appearance	Liquid
Storage	This antibody is stable for at least one (1) year at -20C. Avoid multiple freeze- thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P39748
Dilution Range	Immunoblotting: use at 1-5 ug/ml. Positive controls: MCF-7 cells and recombinant fusion protein.
Expiration Date	12 months from date of receipt.