

Product Datasheet

Anti-cdc34 Antibody (orb1148675)

Description	Anti-cdc34 Mouse Monoclonal Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	N-terminus 6x-His tagged fusion protein encoding full-length human Cdc34 expressed in E. coli. Host Species: Mouse Stabilizers: None Preservatives: None. Available on request.
Target	cdc34
Preservatives	10 mM PBS, pH 7.4.
Form/Appearance	Liquid
Storage	These antibodies are stable for at least one (1) year at -20C to -70C. Avoid multiple freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P49427
Dilution Range	Immunoblotting and Immuno-precipitation: use at 1-5 ug/ml. Positive controls: 293 cells and recombinant fusion protein.
Expiration Date	12 months from date of receipt.