

Product Datasheet

Anti-Nitric Oxide Synthase, Neuronal Antibody (orb1148680)

Description	Anti-Neuronal NO Synthase Rabbit Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Peptide corresponding to aa 19-31 of human neuronal nitric oxide synthase. Host Species: Rabbit Stabilizers: None Preservatives: 0.02% sodium azide
Target	Nitric Oxide Synthase, Neuronal
Preservatives	PBS, pH 7.4.
Form/Appearance	Liquid
Storage	This antibody is stable for 1 month at 4°C or one year at -20°C. Avoid multiple freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Uniprot ID	P29475
Dilution Range	Immunoblotting: use at 1:500-1:1,000 dilution. In immunoblots, a band of approximately 155 kD is detected.
Expiration Date	12 months from date of receipt.