

Product Datasheet

Anti-Protein Kinase D2 Antibody (orb1148681)

Description	Anti-Protein Kinase D2 Rabbit Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Peptide corresponding to aa 849-864 of human PKD2
Target	Protein Kinase D2
Preservatives	0.02% Sodium Azide. PBS, pH 7.4, 0.02% sodium azide.
Form/Appearance	Liquid
Storage	This antibody is stable for at least one (1) year at -20C. Avoid multiple freeze- thaw cycles.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Uniprot ID	Q9BZL6
Dilution Range	Immunoblotting: use at 1:500-1:1,000 dilution. A band of ~105kDa is detected. These are recommended concentrations. Endusers should determine optimal concentrations for their applications.
Expiration Date	12 months from date of receipt.