

Product Datasheet

Anti-NAK/TBK1 (CT) Antibody (orb1148682)

Description	Anti-NAK/TBK1 (CT) Rabbit Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Peptide corresponding to aa 712-727 of human NAK/TBK1. This sequence is identical to that of mouse NAK/TBK1. Host Species: Rabbit Stabilizers: None Preservatives: 0.02% sodium azide
Target	NAK/TBK1 (CT)
Preservatives	PBS, pH 7.4.
Form/Appearance	Liquid
Storage	This antibody is stable for at least one (1) year at -20C. Avoid multiple freeze- thaw cycles.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Uniprot ID	Q9UHD2
Dilution Range	Immunoblotting: use at 0.5-1 ug/ml. In immunoblots, a band of 84 kD is detected. Positive control: MOLT4 cell lysate.
Expiration Date	12 months from date of receipt.