

Product Datasheet

Anti-CIKS/Act1 (CT) Antibody (orb1148684)

Description	Anti-CIKS/Act1 (CT) Rabbit Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Peptide corresponding to aa 554-568 of human CIKS. Host Species: Rabbit Stabilizers: None Preservatives: 0.02% sodium azide
Target	CIKS/Act1 (CT)
Preservatives	PBS, pH 7.4.
Form/Appearance	Liquid
Storage	This antibody is stable for at least one (1) year at -20C. Avoid multiple freeze- thaw cycles.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Uniprot ID	O43734
Dilution Range	Immunoblotting: use at 0.5-1 ug/ml. In immunoblots, a band of 63 kD is detected. Positive control: Human placenta tissue lysate.
Expiration Date	12 months from date of receipt.