

Product Datasheet

Anti-NPFF1 Antibody (orb1148686)

Description	Anti-NPFF1 Rabbit Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	N-terminal sequence MEGEPSQPPNSSWPLS and C-terminal sequence CSHLPLTIPAWDI of human NPFF1. Host Species: Rabbit Preservatives: 0.02% sodium azide.
Target	NPFF1
Preservatives	PBS.
Form/Appearance	Liquid
Storage	These antibodies are stable for at least one (1) year at -20C to -70C. Store product in appropriate aliquots to avoid multiple freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Uniprot ID	Q9GZQ6
Dilution Range	Immunoblotting: use at 1:500-1:1,000 dilution. ELISA: use at 1:1,000-1:5,000 dilution.
Expiration Date	12 months from date of receipt.