

Product Datasheet

Anti-Corticotropin Releasing Factor Antibody (orb1148687)

Description

Anti-Corticotropin Releasing Factor

Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Purified human/rat corticotropin-releasing factor (CRF, Sigma). Host Species: Mouse Antibody Class: IgG1 Purification method: Three sequential precipitations with 40% ammonium sulfate. Purified antibodies yield 1 band in native SDS-PAGE. Preservatives: None. Available on request.
Target	Corticotropin Releasing Factor
Preservatives	PBS, pH 7.4.
Form/Appearance	Liquid
Storage	These antibodies are stable for at least one (1) year at -20C to -70C. Store product in appropriate aliquots to avoid multiple freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P06850
Dilution Range	These antibodies may be used in immunoassays to detect and quantitate CRF. Other applications are under investigation.
Expiration Date	12 months from date of receipt.