

Product Datasheet

Anti-Insulin receptor (beta subunit) Antibody (orb1148696)

Description	Anti-Insulin receptor Mouse Monoclonal Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Human placental insulin receptor
Target	Insulin receptor (beta subunit)
Preservatives	PBS, pH 7.4.
Form/Appearance	Liquid
Storage	This antibody is stable for 1 month at 4°C or one year at -20°C. Avoid multiple freeze-thaw cycles.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P06213
Expiration Date	12 months from date of receipt.