

Product Datasheet

Anti-c-Jun (Phospho-Thr239) Antibody (orb1148717)

Description

Anti-c-Jun (Phospho-Thr239) Rabbit

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Peptide sequence that includes phosphorylation site of threonine 239 (G-E-T(p)-P-P) derived from human c-Jun and conjugated to KLH. Accession no.: P05412, NP_002219.1 Host Species: Rabbit Specificity: This antibody detects endogenous human, mouse, and rat c-Jun only when phosphorylated at threonine 239.
Target	c-Jun (Phospho-Thr239)
Preservatives	0.02% Sodium Azide. PBS (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1.0 mg/mL
Storage	This antibody is stable for at least one (1) year at -20C. Can be stored at 4C for short-term use.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Dilution Range	Immunoblotting: use at dilution of 1:500-1:1,000. A band of ~43kDa is detected. Immunohistochemistry: use at dilution of 1:50- 1:100. These are recommended working dilutions. Enduser should determine optimal dilutions for their applications.
Expiration Date	12 months from date of receipt.