

Product Datasheet

Anti-BAFF-R/BR3 (CT) Antibody (orb1148755)

Description	Anti-BAFF-R (CT) BR3 Rabbit Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Peptide corresponding to aa 169-183 at the C-terminus of human BAFF-R. Host Species: Rabbit Stabilizers: None Preservatives: 0.02% sodium azide.
Target	BAFF-R/BR3 (CT)
Preservatives	PBS, pH 7.4.
Form/Appearance	Liquid
Storage	This antibody is stable for at least one (1) year at -20C. Avoid multiple freeze- thaw cycles.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Uniprot ID	Q5H8V1
Dilution Range	Immunoblotting: use at 2-5 ug/ml. Positive control: human, mouse, or rat spleen tissue lysates. Immunohistochemistry: use at 5 ug/ml. Positive control: human tonsil.
Expiration Date	12 months from date of receipt.