

Product Datasheet

Anti-Secretogranin II Antibody (orb1148860)

Description	Anti-Secretogranin II Rabbit Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Monkey, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Synthetic peptide derived from the N-terminus of mature SgII protein: QRNQLLGKEPDLRLNVQK (Accession no. AAA36607). Host Species: Rabbit Specificity: This antibody recognizes human, rat, and squirrel monkey SgII.
Target	Secretogranin II
Form/Appearance	Liquid
Storage	This antibody is stable for at least one (1) year at -20C to -70C. Store antibody in appropriate aliquots to avoid multiple freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Uniprot ID	P13521
Dilution Range	Immunoblotting: use at 0.5-1.0ug/ml. Immunohistochemistry: use at 2-5ug/ml.
Expiration Date	12 months from date of receipt.