

Product Datasheet

Anti-Caspase-3 Antibody (orb1148884)

Description

Anti-Caspase-3 Antibody Rabbit Polyclonal

Species/Host	Rabbit
Reactivity	Bovine, Canine, Gallus, Guinea pig, Hamster, Human, Monkey, Mouse, Porcine, Rabbit, Rat, Sheep, Yeast
Conjugation	Unconjugated
Immunogen	Recombinant full length human caspase-3. Accession no.: P42574 Gene ID: 836 Host Species: Rabbit Antibody Class: Polyclonal Specificity: This antibody recognizes caspase-3 in human, mouse, rat, bovine, canine, chicken, guinea pig, hamster, monkey, pig, rabbit, and yeast.
Target	Caspase-3
Preservatives	0.09% Sodium Azide. PBS, 50% glycerol, 0.09% sodium azide.
Form/Appearance	Liquid
Storage	This antibody is stable for 1 month at 4°C or one year at -20°C. Avoid multiple freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Uniprot ID	P42574
Entrez	600636
Dilution Range	Immunoblotting: use at 1-5ug/ml. A band of ~36kDa is detected. Immunohistochemistry: use at 10ug/ml These are recommended concentrations. Endusers should determine optimal concentrations for their applications.
Expiration Date	12 months from date of receipt.