

Product Datasheet

Anti-Angiotensin Converting Enzyme 2 (ACE2) (orb1148959)

Description

Anti-Angiotensin Converting Enzyme

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein corresponding to aa 406-677 of human ACE2 protein.
Target	Angiotensin Converting Enzyme 2 (ACE2)
Preservatives	0.01% Thimerosal. 0.1M Tris, 0.1M Glycine, 20% glycerol, 0.01% Thimerosal.
Form/Appearance	Liquid
Concentration	0.73 mg/ml
Storage	This antibody is stable for 1 month at 4°C or one year at -20°C. Avoid multiple freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Uniprot ID	Q9BYF1
Entrez	300335
Dilution Range	Immunoblotting: use at dilution of 1:1,000. A band of ~92kDa is detected. Detection of ACE2 in Hela cell lysate (30ug) with #19505 diluted 1:1,000.
Expiration Date	12 months from date of receipt.