

Product Datasheet

Anti-Angiotensin Converting Enzyme 2 (ACE2) (orb1148960)

Description

Anti-Angiotensin Converting Enzyme

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Synthetic peptide corresponding to amino acids near the N-terminus of human ACE2 protein.
Target	Angiotensin Converting Enzyme 2 (ACE2)
Preservatives	0.02% Sodium Azide. PBS, 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1.0 mg/mL
Storage	This antibody is stable for 1 month at 4°C or one year at -20°C. Avoid multiple freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Uniprot ID	Q9BYF1
Entrez	300335
Dilution Range	Immunoblotting: use at 1ug/ml. A band of ~92kDa is detected. Detection of ACE2 in human kidney lysate with #19504 at 1ug/ml.
Expiration Date	12 months from date of receipt.