

Product Datasheet

Anti-FKBP52 Antibody (orb1149088)

Description

Anti-FKBP52 Mouse Monoclonal Antibody

Species/Host	Mouse
Reactivity	Canine, Hamster, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	Recombinant full-length FKBP52 with an N-terminal polyhistidine tag. Host Species: Mouse Antibody Class: IgG1 Preservatives: 0.09% sodium azide Other Additives: 50% glycerol
Target	FKBP52
Preservatives	PBS, pH 7.4.
Form/Appearance	Liquid
Storage	This antibody is stable for at least one (1) year at -20C. Avoid repeated freezing and thawing.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Uniprot ID	Q02790
NCBI	NP_002005.1
Entrez	2288
Dilution Range	Immunoblotting: use at 0.5-1ug/ml. A band of ~52 kDa is detected. Immunohistochemistry: use at 4-5ug/ml Immunoprecipitation: 5ug on 10-20ul Protein A-Sepharose + 100ul sample. These are recommended concentrations. User should determine optimal concentrations for their application. Positive control: HeLa cell lysate, rat kidney and rat/mouse testes tissue extracts.
Expiration Date	12 months from date of receipt.