

Product Datasheet

Anti-Estrogen Receptor beta Antibody (orb1149140)

Description

Anti-Estrogen Receptor beta Mouse

Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	6-His fusion protein containing the region encoding aa 1-153 of human estrogen receptor- β (ER- β) expressed in E. coli. Host Species: Mouse Antibody Class: IgG2b/ μ ° Stabilizers: None Preservatives: None. Available on request.
Target	Estrogen Receptor beta
Preservatives	PBS, pH 7.4
Form/Appearance	Liquid
Storage	This antibody is stable for at least one (1) year at -20°C . Avoid multiple freeze-thaw cycles.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	Purified
Uniprot ID	Q92731
NCBI	NP_001035365
Entrez	2100
Dilution Range	Immunoblotting: use at 1-5 ug/ml. Positive controls: Native receptor from human thymus, spleen, ovary, and testis or recombinant protein.
Expiration Date	12 months from date of receipt.