

Product Datasheet

Anti-BRCA1 Antibody (orb1149149)

Description

Anti-BRCA1 Mouse Monoclonal Antibody

Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	GST fusion protein expressed in E. coli corresponding to aa 341-748 of full-length BRCA1. Host Species: Mouse Antibody Class: IgG1/ μ Stabilizers: None Preservatives: None. Available on request.
Target	BRCA1
Preservatives	PBS, pH 7.4
Form/Appearance	Liquid
Storage	This antibody is stable for at least one (1) year at -70C. Avoid multiple freeze- thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P38398
NCBI	NP_009225
Entrez	672
Dilution Range	Immunoblotting: use at 1-10 ug/ml. Immunoprecipitation: use at 1-10 ug/ml. Positive controls: Any normal human tissue, especially those containing rapidly proliferating cells, such as lymphoid germinal centers.
Expiration Date	12 months from date of receipt.