

Product Datasheet

Anti-Cyclin E1 Antibody (orb1149150)

Description	Anti-Cyclin E1 Rabbit Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Peptide corresponding to amino acids 43-52 of cyclin E1. Host Species: Rabbit Stabilizers: None Preservatives: 0.02% sodium azide.
Target	Cyclin E1
Preservatives	PBS, pH 7.4
Form/Appearance	Liquid
Storage	This antibody is stable for at least one (1) year at -70C. Avoid multiple freeze- thaw cycles.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Uniprot ID	P24864
NCBI	NP_001229
Entrez	898
Dilution Range	Immunoblotting: use at 1:500-1:1,000.
Expiration Date	12 months from date of receipt.