

Product Datasheet

Anti-Cyclin E2 Antibody (orb1149151)

Description	Anti-Cyclin E2 Rabbit Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Peptide corresponding to amino acids 2-15 of cyclin E2. Host Species: Rabbit Stabilizers: None Preservatives: 0.02% sodium azide.
Target	Cyclin E2
Preservatives	PBS, pH 7.4
Form/Appearance	Liquid
Storage	This antibody is stable for 1 month at 4°C or one year at -20°C. Avoid multiple freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Uniprot ID	O96020
NCBI	NP_477097
Entrez	9134
Dilution Range	Immunoblotting: use at 1:500-1:1,000.
Expiration Date	12 months from date of receipt.