

Product Datasheet

**Anti-Retinoblastoma protein phosphoserine
249/phosphothreonine 252 Antibody (orb1149158)**

Description

Anti-Retinoblastoma phosphoserine

	Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Phosphopeptide region of human Rb protein containing serine 249 and threonine 252.
Target	Retinoblastoma protein phosphoserine 249/phosphothreonine 252
Preservatives	0.05% Sodium Azide. PBS, pH 7.4, 0.05% sodium azide, 1mg/ml rAlbumin.
Form/Appearance	Liquid
Concentration	1.0 mg/mL
Storage	This antibody is stable for 1 month at 4°C or one year at -20°C. Avoid multiple freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Uniprot ID	P06400
NCBI	NP_000312
Entrez	5925
Dilution Range	Immunoblotting: use at 0.25-1ug/ml. Endusers should determine optimal concentration for their application.
Expiration Date	12 months from date of receipt.