

Product Datasheet

Anti-ZIP Kinase Antibody (orb1149168)

Description

Anti-ZIP Kinase Rabbit Polyclonal Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC, WB

Applications
Immunogen

Peptide corresponding to aa 279- 298 of human ZIP kinase (accession no. BAA81746). Host Species: Rabbit Stabilizers: None Preservatives: 0.02% sodium azide.

Target

ZIP Kinase

Preservatives

PBS, pH 7.4.

Form/Appearance

Liquid

Storage

This antibody is stable for at least one (1) year at -20C. Avoid multiple freeze-thaw cycles.

Note

For research use only

Clonality

Polyclonal

Purity

Purified

Uniprot ID
O43293
NCBI
NP_001339
Entrez
1613
Dilution Range

Immunoblotting: use at 1ug/ml. Positive control: HeLa cell lysate.
Immunocytochemistry: use at 10ug/ml. These are recommended concentrations. Enduser should determine optimal concentrations for their applications.

Expiration Date

12 months from date of receipt.