

Product Datasheet

Anti-Caspase-9 (I16) Antibody (orb1149181)

Description

Anti-Caspase-9 (116) Rabbit Polyclonal

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Peptide corresponding to aa 41-56 of human caspase-9 (accession no. P55211). Host Species: Rabbit Stabilizers: None Preservatives: 0.02% sodium azide.
Target	Caspase-9 (I16)
Preservatives	PBS, pH 7.4.
Form/Appearance	Liquid
Storage	This antibody is stable for at least one (1) year at -20C. Avoid multiple freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Uniprot ID	P55211
NCBI	NP_001220
Entrez	842
Dilution Range	Immunoblotting: use at 1ug/ml. These are recommended concentrations. Enduser should determine optimal concentrations for their application. Positive control: Whole cell lysate from HeLa cells.
Expiration Date	12 months from date of receipt.