

Product Datasheet

Anti-Daxx (CT) Antibody (orb1149187)

Description

Anti-Daxx (CT) Rabbit Polyclonal Antibody

Species/Host

Rabbit

Reactivity

Human, Monkey, Mouse

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

Peptide corresponding to aa 722- 740 at the C-terminus of human Daxx (accession no. CAG33366). Host Species: Rabbit Stabilizers: None Preservatives: 0.02% sodium azide.

Target

Daxx (CT)

Preservatives

PBS, pH 7.4.

Form/Appearance

Liquid

Storage

This antibody is stable for at least one (1) year at -20C. Avoid multiple freeze-thaw cycles.

Note

For research use only

Clonality

Polyclonal

Purity

Purified

Uniprot ID
Q9UER7
NCBI
NP_001341
Entrez
1616
Dilution Range

Immunoblotting: use at 1:500-1:1,000 dilution. Immunocytochemistry: use at 10ug/ml. These are recommended concentrations. Enduser should determine optimal concentrations for their applications. Positive control: Whole cell lysate from HeLa cells.

Expiration Date

12 months from date of receipt.