

Product Datasheet

Anti-RICK (CT) Antibody (orb1149230)

Description	Anti-RICK (CT) Rabbit Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Synthetic peptide corresponding to aa 508-522 of human RICK (accession no. AF027706). Host Species: Rabbit Stabilizers: None Preservatives: 0.02% sodium azide.
Target	RICK (CT)
Preservatives	PBS, pH 7.4.
Form/Appearance	Liquid
Storage	This antibody is stable for at least one (1) year at -20C. Avoid multiple freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Uniprot ID	O43353
NCBI	NP_003812
Entrez	8767
Dilution Range	Immunoblotting : use at 1-2ug/ml. Positive control: A431 cell lysate. Immunocytochemistry: use at 5ug/ml.
Expiration Date	12 months from date of receipt.