

Product Datasheet

Anti-CCR5 (NT) Antibody (orb1149276)

Description	Anti-CCR5 (NT) Rabbit Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Peptide corresponding to aa 6-20 of human CCR5. Host Species: Rabbit Stabilizers: None Preservatives: 0.02% sodium azide.
Target	CCR5 (NT)
Preservatives	PBS, pH 7.4.
Form/Appearance	Liquid
Storage	This antibody is stable for at least one (1) year at -20C. Avoid multiple freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Uniprot ID	P51681
NCBI	NP_000570
Entrez	1234
Dilution Range	Immunoblotting: use at 1:1,000-1:2,000 dilution. Positive control: Whole cell lysate from THP-1 cells.
Expiration Date	12 months from date of receipt.