

Product Datasheet

Human IL-17E / IL-25 Protein (orb1149282)

Description

Human IL-17E, His Tag is expressed from

Tyr 33 - Gly 177 (A...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

IL-17E

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 18.8 kDa. The protein migrates as 20-30 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

90%

Source

Human IL-17E, His Tag (orb1149282) is expressed from human 293 cells (HEK293). It contains AA Tyr 33 - Gly 177 (Accession # NP_073626.1).

MW

18.8 kDa

NCBI
NP_073626.1
Expiration Date

6 months from date of receipt.