

Product Datasheet

Human TNFSF11 / RANKL / CD254 Protein (orb1149290)

Description Receptor activator of nuclear factor kappa-B ligand (RANKL), also known as tumor necrosis factor lig...

Reactivity Human

Endotoxins 0.1 EU per µg

Conjugation Unconjugated

Target TNFSF11

Preservatives PBS, pH7.4

Form/Appearance Powder

Storage -20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Tag Native

Note For research use only

Application notes This protein carries no "tag". The protein has a calculated MW of 20.5 kDa. The protein migrates as 28-33 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.


Purity 95%

Source Human TNFSF11 Protein, premium grade (orb1149290) is expressed from human 293 cells (HEK293). It contains AA Gly 64 - Asp 245 (Accession # AAC51762.1).

MW 20.5 kDa

NCBI [AAC51762.1](#)

Expiration Date 6 months from date of receipt.



Immobilized
Human
TNFSF11
Protein,
premi...