

Product Datasheet

Human TLR4 ELISA Kit (orb1149419)

Description	Human TLR4 (Toll-Like Receptor 4) ELISA Kit
Reactivity	Human
Range	0.75ng/ml-24ng/ml
Tested Applications	ELISA
Target	TLR4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	0.1ng/ml
Expiration Date	Please enquire.