

Product Datasheet

Human C1QTNF3 ELISA Kit (orb1149429)

Description	Human C1QTNF3 (C1q and Tumor Necrosis Factor Related Protein 3) ELISA Kit
Reactivity	Human
Range	1.25ng/ml-40ng/ml
Tested Applications	ELISA
Target	C1QTNF3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	0.1ng/ml
Expiration Date	Please enquire.