

Product Datasheet

Human Ang1-7 ELISA Kit (orb1149430)

Description	Human Ang1-7 (Angiotensin1-7) ELISA Kit
Reactivity	Human
Range	5pg/ml-160pg/ml
Tested Applications	ELISA
Target	Ang1-7
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	1.0pg/ml
Expiration Date	Please enquire.