

## Product Datasheet

### Human IL-31 ELISA Kit (orb1149435)

|                            |                                                         |
|----------------------------|---------------------------------------------------------|
| <b>Description</b>         | Human IL-31 (Interleukin 31) ELISA Kit                  |
| <b>Reactivity</b>          | Human                                                   |
| <b>Range</b>               | 1.25pg/ml-40pg/ml                                       |
| <b>Tested Applications</b> | ELISA                                                   |
| <b>Target</b>              | IL-31                                                   |
| <b>Storage</b>             | 2-8°C                                                   |
| <b>Note</b>                | For research use only                                   |
| <b>Application notes</b>   | Sandwich                                                |
| <b>Sample Types</b>        | Cells, serum, tissue, saliva, lavage fluid, urine, etc. |
| <b>Sensitivity</b>         | 0.1pg/ml                                                |
| <b>Expiration Date</b>     | Please enquire.                                         |