

Product Datasheet

Human MUC5AC ELISA Kit (orb1149457)

Description	Human MUC5AC (Mucin 5 Subtype AC) ELISA Kit
Reactivity	Human
Range	2ng/ml-64ng/ml
Tested Applications	ELISA
Target	MUC5AC
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	0.1ng/ml
Expiration Date	Please enquire.