

Product Datasheet

Human CYFRA21-1 ELISA Kit (orb1149474)

Description	Human CYFRA21-1 (Cytokeratin Fragment Antigen 21-1) ELISA Kit
Reactivity	Human
Range	0.25ng/ml-8ng/ml
Tested Applications	ELISA
Target	CYFRA21-1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	0.1ng/ml
Expiration Date	Please enquire.