

Product Datasheet

Human CK-18 ELISA Kit (orb1149475)

Description	Human CK-18/KRT18 (Cytokeratin 18) ELISA Kit
Reactivity	Human
Range	2.5ng/ml-80ng/ml
Tested Applications	ELISA
Target	CK-18/KRT18
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	0.1ng/ml
Expiration Date	Please enquire.