

Product Datasheet

Human MIP-3beta ELISA Kit (orb1149492)

Description	Human MIP-3beta/ELC/CCL19 (Macrophage Inflammatory Protein 3beta) ELISA Kit
Reactivity	Human
Range	15pg/ml-480pg/ml
Tested Applications	ELISA
Target	MIP-3β/ELC/CCL19
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	1.0pg/ml
Expiration Date	Please enquire.