

Product Datasheet

Human THP ELISA Kit (orb1149582)

Description	Human THP (Tamm-Horsfall Glycoprotein) ELISA Kit
Reactivity	Human
Range	6.25ng/ml-200ng/ml
Tested Applications	ELISA
Target	THP
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	1.0ng/ml
Expiration Date	Please enquire.