

Product Datasheet

Human CXCR3 ELISA Kit (orb1149599)

Description	Human CXCR3 (CXC-Chemokine Receptor 3) ELISA Kit
Reactivity	Human
Range	125pg/ml-4000pg/ml
Tested Applications	ELISA
Target	CXCR3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	10pg/ml
Expiration Date	Please enquire.