

Product Datasheet

Human CTX I ELISA Kit (orb1149601)

Description	Human CTX I (Cross Linked C-telopeptide of Type I Collagen) ELISA Kit
Reactivity	Human
Range	0.625ng/ml-20ng/ml
Tested Applications	ELISA
Target	CTXI
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	0.1ng/ml
Expiration Date	Please enquire.