

Product Datasheet

Human ApoC3 ELISA Kit (orb1149618)

Description	Human ApoC3 (Apolipoprotein C3) ELISA Kit
Reactivity	Human
Range	1mg/dl-32mg/dl
Tested Applications	ELISA
Target	ApoC3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	0.1mg/dl
Expiration Date	Please enquire.