

Product Datasheet

Human ApoC1 ELISA Kit (orb1149620)

Description	Human ApoC1 (Apolipoprotein C1) ELISA Kit
Reactivity	Human
Range	0.625mg/dl-20mg/dl
Tested Applications	ELISA
Target	ApoC1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	0.1mg/dl
Expiration Date	Please enquire.