

Product Datasheet

Human ApoA4 ELISA Kit (orb1149622)

Description	Human ApoA4 (Apolipoprotein A4) ELISA Kit
Reactivity	Human
Range	7.5ug/ml-240ug/ml
Tested Applications	ELISA
Target	ApoA4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	1.0ug/ml
Expiration Date	Please enquire.