

Product Datasheet

Human BMPR2 ELISA Kit (orb1149628)

Description	Human BMPR2 (Bone Morphogenetic Protein Receptor II) ELISA Kit
Reactivity	Human
Range	0.625ng/ml-20ng/ml
Tested Applications	ELISA
Target	BMPR2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	0.1ng/ml
Expiration Date	Please enquire.