

Product Datasheet

Human IGFBP-6 ELISA Kit (orb1149634)

Description	Human IGFBP-6 (Insulin-like Growth Factor Binding Protein 6) ELISA Kit
Reactivity	Human
Range	250ng/ml-8000ng/ml
Tested Applications	ELISA
Target	IGFBP-6
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	10ng/ml
Expiration Date	Please enquire.