

Product Datasheet

Human ACE2 ELISA Kit (orb1149652)

Description	Human ACE2 (Angiotensin I Converting Enzyme 2) ELISA Kit
Reactivity	Human
Range	0.3125ng/ml-10ng/ml
Tested Applications	ELISA
Target	ACE2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	0.1ng/ml
Expiration Date	Please enquire.